

AN/PDR-75A Coming Soon

HEY, DID YOU HEAR THE NEW AN/PDR-75A RADIAC SET IS ON THE WAY!

THIS IS REALLY EXCITING, DON'T YOU AGREE?!

YEAH, SURE! I ALWAYS ENJOY BEING REPLACED.

JUST BE AWARE, I HEARD THIS NEW ONE IS A BIT SENSITIVE WHEN IT COMES TO USING THE WRONG WRIST DOSIMETER.

The new AN/PDR-75A radiac set will soon hit the field. It will be a big step up from the old AN/PDR-75 in that it's 1,000 times more sensitive. And it's smaller and lighter.

The -75A will work essentially like the -75. You wear a wrist dosimeter and at the end of operations it's placed in the -75A to measure the radiation dose.

But the -75A's DT-236A/P green wrist dosimeter, NSN 6665-01-602-9991, can be used only with the -75A. And the -75's DT-236/PDR-75 black wrist dosimeter, NSN 6665-01-043-2191, can be used only with the -75. They are **not** interchangeable. If you use the wrong dosimeter, you may get an incorrect radiation reading.

When you get the new AN/PDR-75A, make sure you have enough DT-236A/P wrist dosimeters for your unit.

Refer to TM 11-6665-415-10 and TM 11-6665-415-23&P for more info.

HERE'S HOW TO TELL WHICH DOSIMETER YOU HAVE...

...THE -75A'S IS GREEN, THE -75'S IS BLACK.